
Dräger

Inkubator 8000

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Instructions for use

Inhalt

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Contents

Contents	Page
For Your Safety and That of Your Patients	4
Intended Use	6
Initial preparation	7
Assembly	7
Ventilator Hose Support Arm	7
Mounting Plate for Ancillary Equipment	8
Storage Cupboard	8
Accessories for Rail Mounting	9
Special Accessories	10
Preparing for Operation	11
Air Filter and Water Supply	11
Environmental Sensor Housing	12
Arranging Hoses and Cables	13
Tests of Readiness for Operation	14
Electrical Supply	14
O ₂ Connection	14
Switching On	15
Suppressing Warnings	16
Power Failure Warning and Testing Battery	17
Height Adjustment	18
Visual Inspection	19
Compulsory Calibration	19
Operation	20
Precautions	20
Warming-Up Time	20
Avoiding Overheating	20
Monitoring the Breathing Gas Temperature	22
Avoiding Noise	22
Electrical Safety	22
Fire Risks from Oxygen	22
Physiological Risks from Oxygen	22
Operation	23
Switching On	23
Setting a Comfortable Working Height	24
Placing the Baby in Incubator	24
Adjusting the Baby's Position	24
Note on Setting the Controls	25
Air Temperature Control	26
Skin Temperature Control *)	28
Preparation	28
Operation	29
Control Characteristics	31
Oxygen Supply	31
With Control Module *)	31
Without Control Module	34
Humidity Control	35
Accessories	37
Baby Scales	37
Vacuum Mattress	38
Phototherapy Unit 800	39
Inspection and Maintenance	40
Fault - Cause - Remedy	44

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For Your Safety and that of your Patients

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For Your Safety and that of your Patients¹⁾

For correct and effective use of the apparatus and to avoid hazards it is essential to read the following recommendations and to act accordingly:

Strictly follow the instructions for use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

Maintenance²⁾

The apparatus must be inspected²⁾ and serviced²⁾ by experts at regular 6 month intervals (and a record kept). We recommend obtaining a service contract with DrägerService.

Repairs²⁾ and general overhaul of the apparatus may only be carried out by DrägerService. Only original Dräger spare parts may be used for maintenance. Observe chapter "Maintenance Intervals".

Power connection

The apparatus is to be used only in rooms with mains power supply installations complying with national safety standards (such as in F. R. of Germany: VDE 0107). The requirements laid down in IEC 601/1 "Safety of Medical Electrical Equipment" are applicable for electrically powered equipment.

Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

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Connection with other electrical equipment

Electrical connections to equipment which is not listed in these instructions for use should only be made following consultations with the respective manufacturers or an expert.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevo- cably transferred to the owner or operator to the extent that the apparatus has been serviced or repaired by personnel not employed or authorized by DrägerService or when the apparatus was used in a manner not conforming to its intended use.

Drägerwerk Aktiengesellschaft cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Drägerwerk Aktiengesellschaft are likewise not modified by the recommendations given above.

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1) Insofar as reference is made to laws, regulations or standards, these are based on the legal system of the Federal Republic of Germany.

2) In accordance with DIN 31 051:

Inspection	= examination of actual condition
Service	= measures to maintain specified condition
Repair	= measures to restore specified condition
Maintenance	= inspection, service and, if applicable, repair.

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Intended Use

The 8000 series of incubators are designed for premature and new-born babies.

Controlled temperature, humidity and oxygen concentration*), together with protection from infection, assist even very small babies (<1500 g) to make good progress.

Attention is drawn to possible contra-indications when using skin temperature control*) on babies in shock (p. 28) and during non-physiological use of oxygen for hypoxaemia or hyperoxaemia (p.32).

This equipment is subject to German Medical Regulations and calibration is compulsory before measuring skin temperature (see p. 22).

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
Initial preparation


Assembling the Canopy

After unpacking the incubator, remove all packaging material as well as the transport safety lock on the fan (see "Stripping Down" p.45).

- 1 Fit the canopy at a slight angle, as shown in the illustration.
- 2 Fit the hinged door at an angle and push the leaf spring in the hinge downwards until the door engages in the pivot.

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
Ventilator Hose Support Arm

- Open the hinged door and raise the bed. Insert the cone on the support arm into the hole at the right or left on the far side of the bed and screw tight (using a wide-bladed screwdriver or coin).

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

Mounting Plate for Ancillary Equipment (Dräger-RS-System)

The mounting plate can be attached either on the right or left of the incubator. Maximum load 25 kg

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Fitting Ancillary Equipment to Plate:

- Remove both bottom strips
- 1 Tilt equipment forwards by 45°.
- Insert front latches into the slots in the mounting plate.
- 2 Lower the equipment at the rear. Insert the rear latches into the slots in the mounting plates. Lock with knurled screws.

Storage Cupboards

A cupboard can be attached to the right and/or left of the unit.

- Before fitting the cupboard to the trolley, remove the drawers by lifting them at the front.
- Fit cupboard onto the trolley and secure, following the instructions in the cupboard.
- Replace drawers.
- Fill the drawers with whatever is required.

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
Accessories for Rail Mounting


Accessories such as an O₂ flowmeter, infusion stand, bronchial aspirator or an instrument tray may be attached to the side rail which is also the handle.

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







Preparing for Operation

When all items and components required have been cleaned, disinfected, and fitted as described, proceed as follows:

Air Filter and Water Supply

- 1 Open the rear panel and lower.
- 2 Undo the screws on the filter cover plate.
- 3 Where air filters are already fitted, check the installation date of the filters. If they have been in use for more than two months, replace with new filters. Write the date of fitting on the filters before fitting them.

When fitting: the arrow on the edge of the filter must point towards the incubator. All filter compartments must be filled or filtration will be ineffective.


- 4 Remove empty water bottles from the socket and replace with full bottles.
- Use only infusion bottles which conform to ISO 500 - DIN 58363 - glass 1. These may be refilled with 500 mL clear, distilled or de-mineralised water from the pharmacy. Additives must not be used.

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

Environmental Sensor Housing

(For measuring the temperature, relative humidity and O₂concentration)


- Fit two O₂sensors on incubators with O₂control.
(Use similar expiry dates, if possible)
- Unscrew the retaining ring. Fit O₂ sensor. Replace retaining ring.

- Fit the sensor housing onto the spindle and swing it into place.
- Open the upper flap on the pedestal cover.
- Open the incubator rear panel and lower it.
- Slide the connector from the sensor housing through the slot in the compartment between the filter and bottle compartments. Fit the probe into the socket and secure.

- Close the incubator rear panel.
- Close the pedestal cover.

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


Arranging Cables and Hoses

- 1 Route hose (or cable) through the silicone hose seals.
- 2 Guide hose (or cable) around the support hook on the inside surface.
- 3 Snap the ventilator hoses and cable into the clip on the end of the ventilator hose support arm.

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

Tests of Readiness for Operation

(These instructions apply to incubators fitted with all options. Ignore tests which are not relevant.)

Electrical Supply

- The mains voltage must correspond to the mains voltage specified on the rating plate of the incubator.
- Connect the incubator to the mains supply. Do not push the excess length of cable into the pedestal cover, nor fit it to the upper part of incubator.

1 Lower the flap to see the brief instructions for use.

- Open the doors on the rear panel of the pedestal cover which are held closed by magnets.

2 Make sure that the master switch at the bottom of the support column is switched on (overload circuit breaker).

O₂ connection:


With O₂ control module

3 Use a connecting hose with an M 12 x 1 threaded connector and an oxygen probe. The probe should be in the parking position in the terminal unit.

Without O₂ control module

4 Use a flowmeter with an oxygen probe. The probe should be in the parking position in the terminal unit.

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

Tests of Readiness for Operation

Switching On

- Lower the flap to see the brief instructions for use.
- 1 Switch ON.
- 2 Green LED is lit (LED - Light Emitting Diode).

Whilst the functioning of the microprocessor and its peripherals is being self-tested, dashes will be shown on the actual value displays of the modules for about 35 sec.

If the over-temperature warning is displayed (thermometer symbol) following this self-testing:


Press "reset" button on the skin temperature module.

- 3 Press button - all displays except the power failure LED will light up briefly, a warning tone will sound briefly, and all digital displays will indicate "88.8". All displays will then go blank until the previous settings are shown.

The function test for all lamps, display segments and the audible warning signal can also be carried out during use.


This function test should be carried out at least once a day.

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

Tests of Readiness for Operation

If a fault occurs in one of the control modules "Err" (error) is displayed on the appropriate module. If such a module is switched on, an "inop" warning is also displayed, switching off the entire incubator. Once the faulty module has been switched off, the remaining functions of the incubator will be operational again. Call DrägerService.

- 1 Central warning. The red warning LED flashes at the same time as the individual control warnings.
- 2 Red warning LED "inop" flashes: fault in the central control system of the incubator or in a control module which will also display "Err".

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Suppressing Warnings

- 1 Press button. The acoustic warning will be suppressed for 10 min (does not apply to "inop" and mains failure warnings).
- 2 The yellow pilot LED will light to show that warnings are suppressed.

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

Tests of Readiness for Operation

Power Failure Warning and Testing Battery


- Do not connect the incubator to the mains power supply.
- 1 Press master switch button: An intermittent tone sounds and remains at a constant volume for at least 30 seconds.
 - 2 Red mains failure warning LED is lit.
 - 3 Press the master switch button again. The mains failure warning is cancelled.
- Connect the incubator to the mains power supply.

If the volume of the warning tone drops prematurely:

- Press master switch button to switch on the incubator and leave it switched on for 24 hours to recharge the battery.
- Disconnect from the mains power supply after 24 hours : the volume of the intermittent tone must remain constant for at least 30 seconds.
- Reconnect to mains or press master switch button = intermittent tone ceases.

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
Height Adjustment

Keep the space above the incubator clear of obstructions.


Nothing must be kept on top of the storage cupboard. Hoses and cables must be long enough to allow for the height of the incubator being adjusted, and they must be routed so that kinking, forcible disconnection or squashing cannot occur.

- 1 Press button. The incubator is raised.
- 2 Press button. The incubator is lowered.

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Check that all equipment fitted is ready for operation, in accordance with the relevant "Instructions for Use".

Wear and Tear

The condition of the hinges and pivots of the canopy is very important for the safe operation of the incubator and these should be inspected regularly. If there are any signs of wear and tear, the affected components should be replaced by DrägerService without delay.

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Compulsory Calibration

(FR of Germany only)

The skin temperature sensor and the module for measuring and controlling skin temperature together form an electrical thermometer for determining the temperature of the human body. In West Germany certain legal regulations concerning calibration apply to such thermometer, and users in W. Germany must obtain a copy of the relevant document.

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Operation

(The following instructions apply to incubators fitted with all available options. Disregard the instructions for any options which are not fitted.)

Precautions

Warming-up Time

Allow adequate time for warming-up before placing the baby in the incubator (about 35 min).

Avoiding Overheating

1 Incubator temperature control

The Incubator 8000's temperature control mechanism maintains an extremely stable temperature. Because of the high heating power and good thermal insulation, the inside of the incubator warms up very quickly - even after a malfunction - and for the same reasons cools down slowly when the set temperature is reduced.

Additional external heat sources such as sunshine, heat lamps, radiators, spotlamps and electric blankets (cushions) etc. increase the air temperature inside the incubator in an uncontrolled manner and must be avoided.

2 Incubator air temperature setting

The design and construction of the Incubator 8000 ensure that the baby experiences only minimal heat loss because

- heat loss resulting from radiation is extremely low,
- heat loss resulting from convection is greatly reduced by the low air speed,
- virtually no heat loss occurs from conduction through the mattress and,
- if the relative humidity is set to a high value, heat loss caused by evaporation from the baby's skin is prevented.

By comparison with conventional incubators (e.g. Incubator 6000 or 7000) a lower incubator air temperature can usually be selected with the 8000 series.

The rectal temperature of the baby must be monitored continually - particularly during the first few hours of incubator care.

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3 Reducing the internal temperature of the incubator

To reduce the incubator air temperature, set the temperature on the control module to the lower value which is required. Cooling cannot be accelerated by using a setting lower than the actual value required. The rate of cooling depends on the difference between the temperature of the incubator and the temperature of the room, and, to some extent, on the heating power of the air humidifying system in the incubator. Reducing the room temperature and the relative humidity setting will thus accelerate cooling. In urgent cases the temperature can also be reduced by opening the hinged front door, but if either the front door or the hand ports are open it is absolutely essential that there is continuous supervision so that there is no chance of the baby falling out of the incubator.

4 Phototherapy in the Incubator

Phototherapy lamps are external heat sources supplying heat energy to the baby by absorption of light through the skin. Thus, if the temperature and relative humidity settings in the incubator have not been changed, the baby's body temperature may rise when the phototherapy lamp is switched on. To prevent this happening:

- reduce the incubator air temperature by about 2 °C for 15 minutes before commencing phototherapy;
- set the relative humidity in the incubator to 35 - 45 %;
- reduce the room temperature to at least 50°C below the incubator air temperature. (This temperature is appropriate for the Dräger Heraeus 800 phototherapy unit. Other photo-therapy units, particularly those without a built-in fan, may cause even greater heating of the incubator and therefore the room temperature must then be reduced still further.)

The rectal temperature of the baby must be monitored particularly carefully during phototherapy and the supply of fluids to the baby, e.g. by parenteral infusion, must also be increased in order to compensate for increased water loss.

Do not cover the phototherapy lamp or the incubator canopy with cloths, aluminium foil, or other materials, to boost the phototherapeutic effect. This will cause a build-up of heat since in these circumstances the incubator cannot be cooled by the ambient air, a factor which is absolutely essential for operation.

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Monitoring the Breathing Gas Temperature

The temperature of the breathing gas must be monitored on babies being ventilated since the breathing gas hoses may be warmed by the heated air circulating in the incubator.

Avoiding Noise

Using head boxes or other accessories to deliver pressurised gas may result in unpleasantly high noise levels, as does wear and tear on the bearings of the air circulation motor.

Electrical Safety

Use only ancillary equipment which complies with the regulations of IEC 601/VDE 0750.

Fire Risks from Oxygen

Keep all sources of possible ignition, such as naked lights, matches and cigarettes, well away from the room in which the incubator is in use. Textiles, oils and plastics can very easily catch fire and burn rapidly in an oxygen-enriched atmosphere.

All oxygen valves, fittings and seals must be kept free of oil and grease. Open valves slowly.

Do not use the incubator in the presence of flammable anaesthetic gases or disinfection agents, as there is a risk of explosion.

Do not use any electrical equipment inside the incubator other than equipment or instruments expressly designed for use in incubators.

Flammable liquids, such as alcohol, ether and acetone, must not be used or stored in or near the incubator.

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Physiological Risks from Oxygen

The air in the incubator should only be enriched with oxygen if prescribed by a doctor. It is absolutely essential that such oxygen enrichment is controlled on the basis of the arterially measured oxygen partial pressure in the blood of the baby, as this is the only way of avoiding both hyperoxemia, which might cause damage to the eyes, and hypoxemia which might cause damage to the brain.

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Operation

- Move the incubator to where it is to be used.
- Connect to power supply
- Connect to the gas supply, where required.

If oxygen is required on stand-by only, put the probe in the "parking position" at the terminal.

- 1 Lock the castors in position, if required.
- 2 Lower the flap to see the brief instructions for use.
- 3 Open the hinged front door and lower it.
- Pull the bed right out, if necessary.
- Position a suitable nappy/diaper on the bed and slide it back in. Close the hinged front door.

Switching On

- 1 Switch to ON.
- The green operation LED is lit.



The functioning of the microprocessor and its peripherals is self-tested. Dashes are displayed on the actual value display of the module for about 35 seconds.

- 2 Press button. Function test of all lamps, display segments and the acoustic warning.

All displays except the mains failure LED are briefly lit, a warning tone sounds and all digital displays show "88.8". All displays will then go blank until the measured values and settings are shown.

- Carry out this test **at least once a day**.

The test may be carried out during operation. It does not change or cancel any settings.

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

Operation

Setting a Comfortable Working Height

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
Nothing must be kept on top of the cupboard. Hoses and cables must be long enough to allow for the working height of the incubator to be adjusted and they must be routed so that kinking, forcible disconnection or squashing cannot occur.

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
Placing the Baby in the Incubator

- Put the baby in the pre-heated incubator.
(warming-up time about 35 min.)
- Slide the bed back in position.
- Close the hinged front door and check that the latches are properly engaged.

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Adjusting the Baby's Position

- The bed can be tilted using the handwheels so that the baby's head is either raised or lowered.

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Note on Setting the Controls

The following controls need to be set:

- Air temperature control
- Skin temperature control
- Oxygen concentration control

Adjusting the settings simultaneously adjusts the warning limits. These are as follows:

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
In the event of a warning, **do not alter the settings**, since thus would change the control values; press the warning suppression button, which will suppress the warning for 10 minutes. It should be possible to correct the fault in this time.

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


Operation

Air Temperature Control

Check the baby's rectal temperature at regular intervals. Correct the air temperature setting, if required.

The air temperature control system is automatically activated when the incubator is switched on.

If the skin temperature control system is already switched on:

- 1 Press button = ON
The skin temperature control system will be switched off but it is still possible to measure and display the skin temperature.
- 2 The green operation LED is lit.
- 3 The air temperature inside the incubator is shown alternately with the display "SEt" (request to confirm or reset the value).
 - If the sensor housing has been swung out (as required for O₂ calibration) only dashes will be shown on the display.
- 4 Pre-set value flashes. A pre-set value of 33 °C is shown when switching on the incubator. When switching from skin temperature control to air temperature control the actual air temperature is shown as the set air temperature value.
- 5 Press button briefly:
The set value is confirmed. The displays cease to flash.

Keep button pressed:
The set value is increased in stages (max. 37 °C). Release button when the required set value is shown.
- 6 Press button briefly:
The set value is confirmed. The displays cease to flash.

Keep button pressed:
The set value is reduced in stages (min. 28 °C). Release button when the required set value is shown.
- 7 Green LED:
On: Full heating power
Off: Set temperature exceeded
Flashing: Incubator heating following adjustment of the control system.

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


Operation

Increasing the Temperature Range From 37 ° to 39 °C

The temperature range may only be in-creased above 37 °C if prescribed by a doctor, and if this is done the baby's temperature must be monitored very carefully.

- 1 Press button.
- The yellow pilot LED is lit.
- This alters the warning limit for over-temperature from 38 °C to 40 °C.
- 2 Press button to increase the set val-ue (max. 39 °C)
- 3 Press button to reduce the set value (min. 28 °C)

An increase in the range of set values is cancelled automatically if set values ≤37 °C


- 4 Red warning LED for fan failure.
- 5 Red warning LED for temperature deviations of >±1.5 °C of set value.
- 6 Red warning LED for faulty air tem-perature measurement by environ-mental sensor.
- 7 Red warning LED for overheating.

If this lights:
correct the fault and allow the incubator to cool below the warning threshold.

Then:

- 8 Press button. The incubator should continue to operate normally.

The warning signals are also transmit-ted to the central alarm LED.

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Skin Temperature Control

The skin temperature control system must not be used for babies in shock. The skin temperature of these babies occasionally lies well below the core body temperature, so that using a skin temperature control system could lead to overheating. We would, therefore, advise that for babies in shock the incubator is operated with air temperature control only.

Preparing for Skin Temperature Control

- Establish an appropriate temperature inside the incubator using the air temperature control mode.
- Connect the plug of the sensor to the sensor housing.
- Fit the double-sided adhesive ring to the sensor:
- Remove the backing from the adhesive ring;
- Stick the adhesive ring concentrically on the metal part of the sensor.
- Remove the protective film from the other surface of the adhesive ring.
- Attach the sensor to the skin where the temperature is to be measured with the adhesive ring, and also secure the cable with an adhesive plaster, if necessary.

If the baby is lying on its back, the sensor should be attached to the abdomen between the navel and the sternum. If the baby is lying on its stomach, the sensor should be attached in the region of the kidneys.

The sensor must never be under the baby, since it would then measure core temperature rather than skin temperature and use that as the control value.


Check regularly that the sensor is properly fixed to the baby's skin in case the sensor slips, or becomes detached, and gives a faulty measurement. The air temperature would then increase above the value required by the baby (up to 37 °C).

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



Operation

Using Skin Temperature Control:


Do not switch on the skin temperature control until the sensor attached to the baby has reached skin temperature (about 5 min. or when the display is constant).

- 1 Insert the plug of the sensor into the sensor housing.
- 2 Press button = ON.
The air temperature control is switched off, but the air temperature is still displayed.
(Press button again = to switch back to air temperature control).
- 3 Green LED display is lit.
- 4 The sensor temperature is shown alternately with the display "SEt" (request to confirm or reset the set value).
- When the measured value is outside the range 33 °C - 38 °C the display is as follows:

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

- 1 Set value flashes initially, but then takes over the value of the skin temperature control system (max. 37 °C, min. 35 °C).


- 2 Press button: display must show 36 °C ±0,1 °C.

- 3 Press button briefly:
The set value is confirmed. The display ceases to flash.

Keep button depressed:
The set value is increased in stages (min. 35 °C). Release button when the required set value is reached.

- 4 Press button briefly:
The set value is confirmed. The display ceases to flash.

Keep button depressed:
The set value is reduced in stages (min. 35 °C). Release button when the required set value is reached.

If the skin temperature sensor is going to be temporarily removed, e.g. while the baby is being washed, the air temperature control must first be switched on.


- 5 Red warning LED for set value deviations > ± 0.5 °C.

- 6 Red warning LED for faulty skin temperature sensor.

The warning signals are also transmitted to the central alarm LED.

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
Control Characteristics

The control system is designed to maintain the baby's skin temperature at whatever value has been set. To do this the air temperature in the incubator has to be continually increased or reduced as differences between the measured value and the set value occur. Even when the set value is showing great deviations, the temperature limits of 28 °C or 37 °C are never exceeded.

Practical experience has shown that a baby's skin temperature is subject to sudden changes whilst feeding, sleeping or waking, or being handled. Although adjustments for these deviations are made immediately, as the control is following the set temperature a deviation of a few tenths of a degree is quite normal.

If such small deviations are shown, do not alter the set value.

The set value for skin temperature should only be changed if the rectal temperature shows that correction is necessary.

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

Oxygen Supply

The air in the incubator should only be enriched with oxygen if prescribed by a doctor. It is absolutely essential that oxygen enrichment is controlled on the basis of the arterially measured oxygen partial pressure in the blood of the baby, as this is the only way of avoiding both hyperoxemia, which might cause damage to the eyes, and hypoexemia which might cause damage to the brain.

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
Oxygen Supply with O₂ concentration control:

- Remove probe on O₂ connecting hose from parking position and push into O₂ terminal unit in the pipeline system
- 1 Press button = ON.
(Press again = OFF.)
 - 2 Green operation LED is lit.
 - 3 Yellow pilot LED and "CAL" display flash to request calibration of the sensor.
- Swing the sensor housing out from the incubator. If this is not done within one minute, a sensor warning will be given.

Swing the housing out as far as it will go.

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


Operation

- 1 Press button. Yellow pilot LED is lit.
- 2 During calibration, the display flashes alternately between "CAL" and "- - -". The calibration value of 21% O₂ by vol. (oxygen content of ambient air) is shown after about 45 sec. alternately with "SET" (request to confirm or reset the set value).
- A sensor warning is given if calibration has not commenced one minute after the housing has been swung out.
Press "CAL" button.

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

After calibration:

- Swing the housing back into place. The oxygen concentration inside the incubator is now measured and the actual value shown.
- A sensor alarm is given if the housing has not been swung back into place within one minute after calibration.
Swing the housing back into place.
- 3 The pre-set value (21 Vol. % O₂) flashes.
- 4 Press button briefly:
The set value is confirmed. The displays cease to flash.
Keep button pressed:
The set value is increased in stages (max. 40% O₂ by vol.). Release button when the required set value is reached.
- 5 Press button briefly:
The set value is confirmed. The displays cease to flash.
Keep button pressed:
The set value is reduced in stages (min. 21% O₂ by vol.). Release button when the required set value is reached.

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

Increasing the adjustment range beyond 40 % O₂ by vol.

- 1 Press button. Values up to a maximum of 75% O₂ by vol. can be set.
- 2 Yellow pilot LED is lit.
- Press button to set required value
- 3 The set value is shown.
- If the set value is then reduced again to ≤ 40% O₂ by vol., the extension of the range is automatically cancelled.
- 4 Red warning LED for set value deviation > ± 5% O₂ by vol.
- 5 Red warning LED for oxygen measurement error.

The warning signals are also transmitted to the central alarm LED.

After O₂ therapy:


Remove probe from O₂ terminal unit of pipeline system and replace in parking position.

If any malfunction occurs in the oxygen control module (display "Err"), or there is a sensor warning which cannot be cancelled, and the baby still requires an increased oxygen concentration please refer to the next page.

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
Operation

Oxygen supply without O₂ concentration control module:

- Supply the oxygen from the pipeline system via a flowmeter.
- 1 Connect hose to socket.
- Set the flow rate on the basis of the recommended values in the table below, as required.

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
Monitor the concentration with an additional oxygen analyzer.

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


Humidity Control

Depending on room temperature, high relative humidity ($\geq 70\%$) may lead to condensation inside the incubator canopy. Condensed water will impair the view of the baby and should also be avoided for reasons of hygiene.

- 1 Press button = ON.
(Press again = OFF.)
- 2 Green operation LED is lit.
- 3 Relative humidity inside the incubator is shown, alternately with "SEt" (request to confirm or reset set value).
- 4 Pre-set value (60%) is shown.
- Only dashes will be shown in the centre of the display if the sensor housing is swung out (only necessary for O_2 calibration).
- 5 Press button briefly:
The set value is confirmed. The displays cease to flash.

Keep button pressed:
The set value is increased in stages (max. 85% relative humidity). Release button when the required set value is reached.
- 6 Press button briefly:
The set value is confirmed. The displays cease to flash.


Keep button pressed:
The set value is reduced in stages (min. 35% relative humidity). Release button when the required set value is reached.
- 7 Red warning LED for water shortage. Refill water supply (see p. 36).
- 8 Red warning LED for humidity measurement error.

The warning signals are also transmitted to the central alarm LED.

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
Operation

Refilling the Water Supply

- 1 Open the rear panel and lower it. Wait a few seconds for the water to flow back from the vaporizer unit.

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
- 2 Pull the empty water bottles out of the socket and fit full water bottles. Check that the water bottles are properly secured.


- Only use infusion bottles which conform to ISO 500 - DIN 58363 - glass 1, and refill them with a maximum of 500 mL clear, distilled, and de-mineralised water from hospital pharmacy.
Do not use any additives.

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
**Baby Scales**

(Accessory)

- 1 Adjust the bed to maximum height with both handwheels.
- Open and then lower the front door.

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
- 2 Slide the electronic scales under the bed. Use the ridges on the baseplate to guide the scales.
- Close the front door again.
- Weigh the baby in accordance with the "Instructions for Use" supplied with the scales.
- Remove the scales, if necessary.

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**Vacuum Mattress**

(Accessory)

The vacuum mattress can be formed into any shape, and will keep this shape after it has been evacuated. In this way babies can be supported in extreme positions for special treatment. The standard mattress may also remain inside the incubator.

- Open the front door and then lower it to insert the mattress.
- Position the vacuum mattress and form it into the shape required.
- Position the baby on the mattress and mould the mattress around the baby.
- Evacuate the mattress.
- Close valve and remove hose.
- Close the front door.

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
**Phototherapy Unit 800**

(Accessory)

- The trolley carrying the phototherapy unit must be pushed under the incubator.
- Adjust the height of both the phototherapy unit and the incubator to ensure that the under-side of the phototherapy unit is located directly above the top surface of the incubator canopy.

Be careful not to raise the incubator any further now.

- Operate the phototherapy unit in accordance with the relevant "Instructions for Use".

Note the safety advice in section 4 on page 21 about avoiding increased temperatures.

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O₂ sensor capsule : In the event of a malfunction, always replace both capsules. Ensure than both expiry dates are similar. (See "Fault - Cause - Remedy", p. 44).

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Ambient air filter : Replace after two months (See "Preparing for Operation", p. 11.)

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Hand port seals : The seals must be replaced when the material becomes brittle or tacky, or when a proper seal can no longer be guaranteed.

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Maintenance work should be carried out by trained staff only:

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Fault - Cause - Remedy

Master Switch front panel with height adjustment

Fault	Cause	Remedy	
Symbol	Description		
● Alarm	Central warning: Red LED and intermittent tone.	A module has signalled a fault.	Lower flap to see brief instructions for use and locate fault in appropriate module. Refer to "symbol" column for further details.
● "Inop"	Red LED and continuous tone.	Central fault in the electronic control system.	Switch incubator off and then on again. Renewed "inop" message shows incubator not operational. Call DrägerService! If no "inop" message occurs, set required values. Press "reset" button.
● Powersymbol 	Red LED and continuous tone.	Power failure.	1) Insert mains plug into socket. 2) Switch on circuit breaker in the base. 3) Inform hospital technician. 4) Call DrägerService!
● "Inop" Display "Err"	Red LED and continuous tone. "Err" on actual value display of a module.	Fault in the control module.	Switch off the relevant control module. All other functions should remain operational. Call DrägerService!
Arrow symbols  	No adjustment when the switch is operated.	Thermal cut-out has been activated. The unit is already in final position.	Wait until the motor has cooled down. No further height adjustment is possible.

Air temperature control module

Fault	Cause	Remedy
Symbol	Description	
● "±1,5 °C"	Red LED and intermittent tone.	Temperature deviates from set value by more than ±1.5 °C. Below set value: Close the openings in the canopy. Wait for incubator to warm up. Above set value: Avoid additional heat sources, such as radiators, sunshine and lamps. Wait until incubator has cooled down.
● Thermometer symbol 	Red LED and intermittent tone.	Temperature above 38 °C or above 40 °C *) Avoid additional heat sources, such as radiators, sunshine and lamps. Wait until incubator has cooled down. Press "reset" button when the actual temperature dropped below the warning limit.
● Fan symbol 	Red LED and intermittent tone.	Fan blade not fitted. Fan not rotating or rotating too slowly. Check that fan and drive shaft are firmly coupled.
● "Sensor"	Red LED and intermittent tone.	Sensor housing disconnected, faulty or not swung back into place. Incubator temperature below 5 °C Plug the sensor housing firmly into its socket; replace it or swing it back into the canopy, as necessary. Allow the incubator to warm up thoroughly.

*) For set values up to 37 °C switch-off occurs at 38 °C; and for set values of 37,1 bis 39 °C switch-off occurs at 40 °C.

Skin temperature control module

Fault	Cause	Remedy
Symbol	Description	
● "±0,5 °C"	Red LED and intermittent tone.	Temperature deviates from set value by more than ±0,5 °C. Below set value: Check that the sensor is properly fixed to the baby. Above set value: Measure the core temperature of the baby rectally. Inform the doctor responsible for the baby without delay!
● "Sensor"	Red LED and intermittent tone.	Plug not properly inserted into the sensor housing. Faulty sensor. Check connection and correct if necessary. Replace sensor.
Display "Err"	Continuous tone.	Fault in control module. Switch off control module. All other functions should remain operational. Call DrägerService!
Press "Check" Button : Display NOT 36 °C ± 0,1 °C	Temperature reading outside tolerance level.	Temperature sensor not operating. Call DrägerService! Have equipment recalibrated.

Oxygen control module

Fault		Cause	Remedy
Symbol	Description		
● "±5 Vol. %"	Red LED and intermittent tone	Oxygen concentration deviates from set value by more than $\pm 5\%$ by vol.	Below set value: Close openings on canopy. Check oxygen pipeline to piped medical gas system. Above set value: Switch off other O ₂ sources. If O ₂ concentration does not drop disconnect from piped medical gas system and continue operating with manual control. Call DrägerService!
● "Sensor"	Red LED and intermittent tone	Sensor not connected.	Connect sensor.
		Sensor swung out.	Swing sensor back into place properly.
		Sensor housing incorrectly connected.	Check connections and correct, if necessary.
		Incorrect sensor housing. (Without O ₂ -measuring)	Replace incorrect sensor housing.
		Faulty sensors	Replace faulty sensors and then calibrate sensors.
● "Sensor"	Red LED and intermittent tone. "CAL" flashes.	Calibration not started.	Calibrate sensors (p. 31).
● "Cal" Display "CAL"	Yellow LED and intermittent tone. Display flashing.	Requesting recalibration of O ₂ sensors during continuous operation every 24 hours and during startup.	Calibrate sensors (p. 31).
Display "Err"	Continuous tone.	Fault in control module	Switch off control module. All other functions should remain operational. Call DrägerService!
● "Inop" O ₂ -Display	Red LED and continuous tone. O ₂ -actual value flashes.	O ₂ increase after O ₂ control module is switched off.	Disconnect both from piped medical gas system and other O ₂ sources. Switch incubator off and then on again. If fault persists call DrägerService!

Humidity control module

Fault		Cause	Remedy
Symbol	Description		
● "H ₂ O"	Red LED.	Water shortage.	Replace empty water with full bottles.
● "Sensor"	Red LED and intermittent tone.	Faulty sensor.	Switch off control module. Call DrägerService!
		Plug not properly inserted into the sensor housing.	Check connections and correct, if necessary.
Display "Err"	Continuous tone.	Fault in the control module.	Switch off control module. All other functions should remain operational. Call DrägerService!

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Cleaning - Disinfecting

The incubator must be thoroughly cleaned and disinfected after each baby has used it and about once a week. For the cleaning and disinfecting of accessories, see the relevant "Instructions for Use".

The incubator must first be stripped down.

Stripping-down


- Switch off incubator. Disconnect from mains and disconnect the oxygen supply.
- Remove all accessories
- 1 Disconnect skin temperature sensor from housing.
- 2 Open the front door.
- 3 Remove silicone hose seals.
- 4 Remove seals from hand ports.
- 5 Remove the mattress from the bed.
- 6 Lift off the bed. Undo the screw on the ventilator hose support arm using a widebladed screwdriver or coin. Remove the support arm.
- 7 Pull the metal brackets supporting the bed up and out.
- 8 Remove the base.

Caution! Be careful not to burn yourself. The heating system will still be very hot after operation. When the incubator is closed it takes about one hour for the temperature to cool to 70 °C.

- 9 Firmly hold fan and undo knurled nut. Remove the fan.
- 10 Open the rear panel.
- 11 Remove the water bottles (see p.11).
- 12 Disconnect the hose. Unscrew the bottle fixture.

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


In addition, after two months" use:

- Swing the sensor housing out from the canopy and remove it.
- 1 Tilt the incubator canopy slightly backwards but not as far as it will go. Grip from the rear and lift canopy off in the direction of the arrow.
 - 2 Remove the trough

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Various plastics such as plexiglass and macrolon are used in the incubator. These plastics are sensitive to high temperatures ($> 70^{\circ}\text{C}$) and to certain cleaning and disinfecting agents (e.g. alcohol), and exposure to them may cause stress cracks in particular. To prevent any damage, we recommend that only the cleaning agents and disinfectants specified are used.

Cleaning the Canopy and Incubator

Wipe all surfaces with a damp cloth which has been soaked in detergent. Lower the inner panels of the canopy and wipe the inside carefully.

Use a disinfectant for cleaning if there is any risk of infection.

Disinfecting the Canopy and Incubator

Wipe all surfaces of the canopy, dismantled parts, and the unit itself using a suitable disinfectant. We recommend the following disinfectants. Please note the specified maximum concentrations:

Bacillotox®	required concentration up to 2 %
Buraton®	as specified by the manufacturer
Dismozon®pur	as specified by the manufacturer
Gevisol®	required concentration up to 5 %

Allow time for disinfection, then wipe with a clean, damp cloth and dry well.

Sterilize the water bottles and bottle fixture at a temperature of 134°C and sterilize the water hose at a temperature of 120°C (glove programme).

The sensor housing must never be treated with spray or disinfected in a bath.

Do not disinfect with UV radiation as this may cause damage to the plexiglass components (cracking). This risk does not occur with long-wave radiation during phototherapy.

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
Re-assembling the Incubator

When re-assembling, follow the stripping-down process in reverse order.

- 1 Replace the trough.
- Fold up the inside panels of the canopy until the click-stop. Grip incubator canopy from the rear and insert at an angle. Lower the canopy. Re-insert the seals on the hand ports (note left and right) and the silicone hose seals. Close the canopy.
- Insert sensor housing into the incubator canopy.
- Grease the motor shaft (Vaseline).
- 2 Fit the fan. Grip firmly and tighten the knurled not.
- Check that the fan is firmly secured by pulling on it.
- 3 Fit the base.
- 4 Fit the metal brackets for the bed from above. Note the preferred position.
- 5 Hook the bed onto the rear of the metal brackets and position.
- 6 Place the mattress on the bed.
- Open the rear panel.
- 7 Fit air filter (replace if necessary).

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Do not fit any new water bottles yet. Do not add sterile water until the system has been tested for operational readiness immediately before a baby is placed in the incubator.



Ensure that the correct types of sensors are fitted (see information on the plate on the top of the sensor housing).

- Finally, depending on the disinfectant used, switch on the incubator and operate for a few hours to eliminate all traces of disinfectant.
- Check that the incubator is ready for operation immediately before placing a baby in it (see p. 14).

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
Front panel


- 1 Humidity control
- 2 Oxygen concentration control
- 3 Skin temperature control
- 4 Air temperature control
- 5 On/Off switch, height adjustment, and self-test

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On/Off panel


- 1 Fault indicators and warning status
- 2 Warning tone suppression
- 3 Display/indicator test
- 4 Height adjustment
- 5 On/Off switch

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
Description

A special air flow system has been developed for the 8000 series of incubators (see diagram above). This system routes heated air between the double walls on the entire front (hinged front door), over the roof of the canopy and then draws it down again between the double walls at the rear by means of an extraction system. Thus the baby lies in an area with very low air flow so that heat loss by convection is reduced.

The all-round double wall design of the canopy also reduces the loss of body heat by radiation. When the hinged front door and the hand port doors are open, there is still an effective hot air curtain which prevents appreciable cooling inside the canopy.

A new type of mattress gives very low heat loss by conduction through the bed, resulting in a high level of comfort for the baby. It is made of a particularly soft foam plastic, smoothly encased in film to provide an absolutely hygienic seal.

The incubator is equipped with easy-moving, swiveling castors which permit easy movement around the hospital. Two of the castors can be locked. There are rails, which can also be used as handles, on the sides for holding ancillary equipment.

The sensor housing which is located near the baby's head responds very rapidly to changes in the environment inside the incubator and contributes to a very sensitive control system. The housing contains the air temperature sensor, the optional oxygen sensor and the optional relative humidity sensor. It also contains the connection point for the skin temperature sensor.

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The optional oxygen control module permits the required oxygen concentration within the canopy to be pre-selected and controlled automatically. If the actual value deviates from the set value, warning signals are given. The control system permits a set oxygen concentration to be achieved within a very short time (40 % in 1.5 minutes and 60 % in ≤ 7 minutes). This permits rapid action in emergencies and rapid return to a selected oxygen concentration after opening the hinged front door or the hand ports.

The optional humidity control permits a required humidity value to be pre-selected and maintained. The water required for this is supplied from sterile water bottles. Evaporation of the water guarantees completely hygienic air humidification.

In accordance with hygienic requirements, all parts of the incubator which come into contact with the gas breathed by the baby can be disinfected thermally and can be removed from the basic unit.

The electrically-operated height adjustment mechanism permits an optimum working height to be set for both tall and short nurses. In addition, the incubator can be lowered so that the nurse can carry out lengthy procedures whilst sitting down, and there is also adequate legroom beneath it.

The hand port openings have been designed to provide better access for the arms of the nursing staff than round openings do, yet have a smaller opening cross-section, thus reducing heat loss.

The incubator is controlled by a microprocessor. The incubator carries out a general self-test when it is switched on and automatically every 10 minutes thereafter. This test checks the various assemblies in the incubator which are crucial to safety by operating them briefly. The operator is automatically requested to recalibrate the oxygen sensors after switching on the O₂ control system and, thereafter, after every 24 hours of operation. Measured value deviations or system failures are indicated visually and audibly. The heating system and valves are automatically switched off if operating conditions are outside permitted limits.

The control systems for air temperature, skin temperature and humidity are PI-controls (proportional integral controls) and give optimum control.

An additional built-in air-cooler produces rapid cooling as soon as the actual value of the air temperature is a few tenths of a degree above the set value.

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Safety Features

A self-test is conducted when the incubator is switched on. This involves a test of all memories in the microprocessor control system and a check to establish whether the various programme segments are running correctly.

The function of actuators, acknowledgement signals, displays and indicators is checked by switching them on **and** off. For example, the functioning of the fan is checked not by establishing whether power is reaching the electric motor, but rather by measuring the speed of rotation directly on the fan. The wiring to peripheral devices is provided in duplicate.

This complete test is also repeated every 10 minutes during operation.

All options, whether they are activated or not, are tested. An error message is thus indicated for a faulty module even if it is switched off. This ensures that the operator is always fully aware of the actual condition of the incubator.

When the incubator is switched on again, or when a module is switched on, the relevant required value is set to "starting value". The set value display then flashes requesting confirmation or modification of this value. The actual value indication is displayed alternately with the prompt "SEt".

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Warnings

The incubator has a graduated system of warning signals. Faults which occur are signalled in order of importance. If a non-essential option fails, it does not mean that the incubator has to discontinue any functions which are of vital importance.

The set value deviation and water shortage warnings are signalled by an intermittent warning tone which can be suppressed for 10 minutes. Faults which entail greater potential danger are signalled by a continuous tone which cannot be silenced. In addition, the relevant warning LEDs flash to enable the operator to locate the fault rapidly. Each individual warning also activates an indicator on the central warning LED so that a visual signal can still be seen even with the flap showing the brief instructions for use closed up. The warning signals are described in detail below:

1 Warning for set value deviations and water shortage:

Deviations outside the set limits for oxygen concentration ($> \pm 5\%$ by vol. O₂), skin temperature ($> \pm 0,5$ °C) and air temperature ($> \pm 1,5$ °C), in addition to over-temperature (> 38 or 40 °C) and water shortage result in the following warning signals:

- o The actual value display of the relevant control module flashes.
- o The specific warning LED flashes.
- o The central warning LED flashes.
- o Intermittent warning tone sounds (can be suppressed for 10 min.).

The warning tone for set value difference - 1.5 °C air temperature - is suppressed automatically for 30 minutes when the incubator is switched on to avoid a continuous warning signal sounding during the warming-up phase.

The "Reset" warning button for over-temperature has to be reset by the operator to ensure full awareness of the error. All other alarms are reset automatically when the relevant conditions have returned to normal. If a warning has been suppressed, the red alarm LED, which had previously flashed, stays alight but is constant, and a yellow LED indicator above the switch-off button indicates the switched-off mode.

If another warning occurs when warnings are suppressed, the audible warning is re-activated immediately. Depending on the type of fault, this later warning can be suppressed again. The period for automatic reactivation of the first warning tone is then prolonged by the time interval between the two warnings.

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2 Sensor warnings, fan failure warning:

Faults in the sensors for oxygen, skin temperature, humidity, air temperature and failure of the fan produce the following alarm signals:

- Three dashes ("---") flash on the relevant actual value display.
- The actual air temperature value flashes in the event of fan failure.
- The specific warning LED flashes.
- The central warning LED flashes.
- Continuous warning tone (non-suppressible)

When there have been sensor warnings for oxygen concentration or relative humidity and the relevant module has been switched off, the incubator can continue to operate with the other functions. The same thing applies when there is a skin temperature control warning, though, in this case, the faulty sensor must be removed.

3 Module faults:

Warnings for faults in the control modules for oxygen concentration, skin temperature and relative humidity produce the following warning signals:

- The letters "Err" (= Error) appear on the relevant actual value display, flashing if the module is switched on and not flashing if the module is switched off.
- The red "inop" LED is lit when the module is switched on.
- Continuous warning tone (non-suppressible) if the module is switched on.

The functions of the other modules which are operating correctly are retained. The operator can carry on using the incubator after the faulty module has been switched off. Call DrägerService.

4 Total shutdown

Faults in the central microprocessor control system and in the air temperature control module result in complete shutdown of the incubator. The following warning signals are then given:

- Red "inop" LED is lit.
- Continuous warning tone (non-suppressible).

Individual checks should be carried out to determine whether the fault involves temperature. If the incubator shows the "inop" signal after it is switched on again, it is inoperative. It must not be used and must be clearly marked faulty. Call DrägerService.

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Ambient humidity : 15 to 95% relative humidity
(no condensation)

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Warming-up time : 35 minutes from 20 °C to 31 °C with or without air humidification

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Humidification : Evaporation of distilled water by boiling

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Air velocity : About 8 cm/sec over the bed

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Fresh air supply rate : 30 L/min
 Tilting angle of bed : Right and left continuously adjustable up to 10° tilting angle.

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Noise level : < 50 dB (A) inside the canopy.
 Fresh air filter : Particle class P2 in accordance with DIN 3181 factor 2 %

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Classification :  Type B
 The incubator complies with DIN-VDE 0750, parts 1 and 212

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RFI suppression : The 8000 series of incubators are RFI-suppressed in accordance with relevant German regulations. The appropriate German PTT authorities have been informed that this incubator has been marketed and the German PTT authorities are entitled to inspect incubators in the series to establish whether the regulations are observed.

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Identification : The serial no., the part no. of the incubator and the registration no. are given on the rating plate which is fixed to the rear of the trolley next to the power cable entry point. The scope of the equipment is defined by the part number. (See Order List, p. 62).

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

Skin Temperature Control:

Sensor : Only for earthed skin temperature sensors produced by Dragerwerk AG (see Order List, p.39)

	Display section	Skin temperature sensor
Registration symbol:	15.22 87.06	15.21 86.09
Approval certificate no.:	IB.13-78/87	IB.13-61/86
Calibration error limits:	±0.2 °C	±0.1 °C
Calibration error limits of the entire measuring system over the measuring range: ±0.3 °C.		

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These instructions for use apply
only to Incubator 8000 with se-
rial no.:

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Without entry of serial-no. by
Dräger these instructions for
use are provided for general in-
formation only and are not des-
tined to be used with a specific
device!

Drägerwerk
Aktiengesellschaft
Lübeck
Postfach 1339
Moislinger Allee 53-55
2400 Lübeck 1
Tel. (0451) 882-0
Telex 26807-0
Telefax (0451) 882-2080

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© Drägerwerk Aktiengesellschaft